

ALECK INTERNATIONAL




An ISO 13485:2016 Certified & CDSCOGovernmentofIndiaApprovedCompany


AN EXPERT SOLUTIONS FOR



GASTROENTROLOGY

1. MULTIBANDLIGATOR





HEMO CLIP- DISPOSABLE CARTRIDGE



9.

- Cartridge opening 90, 135 degree.
- Tolerance of cartridge opening >4mm.
- Sterile and single use only.

10. JEJUNAL FEEDING TUBE



- Feeding tube size diameter 12, 14, 16 fr.
- Operating length 120 cm, 155 cm.
- Wire guide Accordant 0.038 inch.
- Sterile and Single use only.

11. ESOPHAGEAL DILATOR



- Accessible sizes 5 mm to 20mm outer diameter.
- Sterile / Non sterile.
- Length : 85cm, 100cm.
- Accordant wire guide 0.038 inch.

12. SPRING TIP GUIDE WIRE



- Size Diameter 0.038 inch.
- Length 250 cm.
- Tip : Straight.
- Sterile and Single use only.

13. POLYPECTOMY SNARE

- Loop Diameter 10 to 50 mm.
- Design of Snare Oval, Duckbill, Hexagonal, Crescent.
- Length 200cm.
- Sterile and Single use only.

14. PEG KIT



- Size 20, 24, 26 Fr With Pull technique Facility.
- Snare For Pulling Tube Also Accessible with kit.
- Sterile and Single use only.

15. GASTROSTOMY REPLACEMENT TUBE



- Size : 24 Fr.
- Sterile and Single use only.

16. ENDO CLAMP DEVICE




- Operating length 230 cm.
- Sheath diameter 2.5 mm.
- Clamp opening diameter 30mm / 40 mm &
- Re-Positionable.
- Sterile and Single use only.

17. MEMORY BASKET




- Design Of Basket Oval, Helical, Hexagonal & Tipless.
- Retain Shape for Long Lasting Time.
- Wire Guided Memory Basket also Accessible.
- Device not recommended for lithotripsy.
- Sterile and Single use only.



- Specially Devloped for enhanced protaction and safety
- Design to protect against mucosal damage during retrieval of sharp Foreign body object from esophagus.

E.R.C.P

19. SPHINCTEROTOME




NEEDLE KNIFE

20.

- Three lumen with tapered tip.
- 200 cm Operating length & total length 230 cm.
- Accordant with 0.035 inch wire guide.
- Cutting Loop with protective insulation sheath.
- Cutting Loop 20/25/30 mm.
- Monopolar cable Accessible Separately.
- Length ,Cutting loop and tip are Accessible on request.
- Design of Knife Tapered Tip, Hook Type, Ball Tip Knife
- Length of cutting wire 7 mm.
- Working length 200 cm ; Total Length 230 cm.
- Sterile and Single use only
- 21. **NEEDLE KNIFE WIRE GUIDED**



- Three Lumen with tapered tip.
- Length of cutting wire 7 mm.
- Working length 200 cm ; Total Length 230 cm.






- Design of Tip Tapered tip, Double lumen, Ball tip
- Tip Diameter 4.5, 5.5 Fr.
- Length 230 cm.
- Wire guided diameter 0.035 inch.
- Sterile and Single use only

26. STONE EXTRACTION BALLOON




- Triple lumen / Double lumen.
- Sheath diameter 7 Fr.
- Accordant with 0.035 inch wire guide.
- Operating length 200 cm ; Total Length 230 cm.
- Sterile and Single use only.

27. STONE RETRIEVAL BASKET



- Double / Single Lumen .
- Basket Diameter 15, 20, 25, 30, 35 mm.
- Length 230 cm.
- Accordant with 3.2 mm.
- Sterile and Single use only.



32. BALLOON INFLATION DEVICE (MICRO)



- Controlled push of air/water.
- 0 to 30 ATM indicator gauge meter.
- 30ml syringe with 3 way stop-cock.
- Reusable and non-strile.
- Used for Esophageal & Bliary Balloon Dilator.

33. BALLOON INFLATION DEVICE (PNEMATIC)



- Reusable, non-sterile.
- For Achalasia Balloon Dilatation.


34. CYTOLOGY BRUSH WIRE GUIDED



35. CYSTOSTOME



- Single Lumen / Double Lumen.
- Accordant with 3.2 mm.
- Length 170/ 230 cm.
- Wire Guide accordant 0.035
- Sterile and single use only.
- Size : 6, 8.5, 10 Fr.
- Length 200 cm.
- Accordant with 2.8 mm.
- Sterile and single use only.



40. BALLOON DILATOR- 3 STEP



- Expanded balloon diameter
- Size 10-11-12 / 12-13-15 / 15-16-18 / 18-19-20
- Balloon length : 5cm.
- Operating length : 180cm , 240cm.
- Wire guide Accordant 0.035 inch.
- Sterile and Single use only

41. BILIARY BALLOON DILATOR



- Expanded balloon diameter 6mm, 8mm, 10mm.
- Balloon length : 4cm.
- Shaft diameter : 2.5mm.
- Operating length : 200 cm.
- Wire guide Accordant : 0.035 inch.
- Sterile and Single use only.



- Accordant with 2.8 mm instrument channel.
- Operating length : 180 cm.
- Atraumatic Prong to reteriave deflated Esophageal Gastric Balloon.
- Sterile and Single use only.

43. GASTRIC BALLOON DEFLATION DEVICE



- Operating Length 170 cm.
- Sheath Diameter 2.5 mm.
- Sterile and Single use only.

BRONCOSCOPY

44. CHENNEL CLEANING BRUSH



- Reuseble / Disposable .
- Langth 120 cm, 200 cm.
- Brush Length : 1.5 cm.
- Accordant with 2.0, 2.8 mm instrument channel.
- Tip : Atraumatic.
- Non Sterile.



- Accessible in 19, 20, 21 swg.
- Operating length : 120 cm.
- Sterile and Single use only.



- Single Lumen.
- Accordant with 2.0, 2.8 mm.
- Length 120, 170 cm.
- Sterile and Single use only
- 47. ACCESSORIES HANGER
- 48. MOUTH GUARD DISPOSABLE WITH OR WITHOUT STRAP
- 49. ENDOSCOPE HANGER
- 50. **BIOPSY VALVE**
 - 51. COMBO HYDRYO BLACK COBRA WIRE

Contact us

ALECK INTERNATIONAL

A-510 Tirupati Plaza, Near Collector Office,

Athwagate,Surat-395001,Gujarat-INDIA

Mobile No: +91 8200537380 aleckinternational12@gmail.com

www.aleckinternational.com