

# STAPLING SOLUTIONS

## Product Reference

### Signia™ Stapling System

MEDTRONIC

Length	Device	Reorder Code	# Procedures
—		SIGPHANDLE	300
Short		SIGADAPTSHORT	50
Standard		SIGADAPTSTND	50
Long		SIGADAPTXL	50

ETHICON

EES Length	EES Product Code
—	—
Short	PCEE45A & PCEE60A
Standard	PSEE45A & PSEE60A
XL	PLEE45A & PLEE60A

### Endo GIA™ Ultra Universal Handle

MEDTRONIC


Length	Device	Reorder Code	Reloads
Short		EGIAUSHORT	25
Standard		EGIAUSTND	25
XL		EGIAUXL	25

ETHICON

EES Length	EES Product Code
Short	SC45A & SC60A
Standard	EC45A & EC60A
XL	EC45AL & Long60A

### Endo GIA™ Reloads with Tri-Staple™ Technology

MEDTRONIC

Standard Reload Length	Device	Reorder Code	Color Code	Staple Heights	Tissue Type	EES Product Code	EES Reload Color
30 mm		SIG30AV†	Gray	2 mm, 2 mm, 2 mm	Vascular		
		SIG30AVM	Tan	2 mm, 2.5 mm, 3 mm	Vascular/Medium		
		SIG30AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick		
45 mm		EGIA45AV†	Gray	2 mm, 2 mm, 2 mm	Vascular	ECR45M	Gray
		EGIA45AVM	Tan	2 mm, 2.5 mm, 3 mm	Vascular/Medium	GST45W	White
		EGIA45AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	GST45B	Blue
		EGIA45AMT	Purple	4 mm, 4.5 mm, 5.0 mm	Medium/Thick	GST45D	Gold
		EGIA45AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	GST45G	Green
		SIG45AXT	Black	4 mm, 4.5 mm, 5 mm	Extra-Thick	GST45T	—Black
60 mm		—	—	—	—	ECR60M	Gray
		EGIA60AVM	Tan	2 mm, 2.5 mm, 3 mm	Vascular/Medium	—	—
		EGIA60AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
		EGIA60AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
		EGIA60AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
		SIG60AXT	Black	4 mm, 4.5 mm, 5 mm	Extra-Thick	—	—
Reinforced Reload Length	Device	Reorder Code	Color Code	Staple Heights	Tissue Type	EES Product Code	EES Reload Color
45 mm		SIGTRSB45AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
		SIGTRSB45AXT	Black	4 mm, 4.5 mm, 5 mm	Extra-Thick	—	—
60 mm		SIGTRSB60AMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
		SIGTRSB60AXT	Black	4 mm, 4.5 mm, 5 mm	Extra-Thick	—	—
Signia™ Small Diameter Reload	Device	Reorder Code	Color Code	Staple Heights	Tissue Type	EES Product Code	EES Reload Color
30 mm		SIGSDS30CTV	Gray	2 mm, 2 mm	Vascular	—	—
		SIGSDS30CTVT	White	2.5 mm, 2.5 mm	Vascular/Medium	VASECR35PEV35A	White
45 mm		SIGSDL45CTVT	White	2.5 mm, 2.5 mm	Vascular/Medium	—	—

†While offering a slimmer fixed anvil, reloads is a vascular reload that does not come with Tri-Staple™ technology.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.  
 †\*\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.  
 07/2021-US-ST-2100120-[WF#2835774]


# STAPLING SOLUTIONS

## Product Reference


### Endo GIA™ Reloads with Tri-Staple™ Technology

MEDTRONIC						ETHICON	
Curved Tip Reload Length	Device	Reorder Code	Color Code	Staple Heights	Tissue Type	EES Product Code	EES Reload Color
30 mm		SIG30CTAV†	Gray	2 mm, 2 mm, 2 mm	Vascular	—	—
		SIG30CTAVM	Tan	2mm, 2.5mm, 3 mm	Vascular/Medium	VASECR35/PVE35A	White
45 mm		SIG45CTAV†	Gray	2 mm, 2 mm, 2 mm	Vascular	—	—
		SIG45CTAVM	Tan	2 mm, 2.5 mm, 3 mm	Vascular/Medium	—	—
		SIG45CTAMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
60 mm		SIG60CTAVM	Tan	2 mm, 2.5 mm, 3 mm	Vascular/Medium	—	—
		SIG60CTAMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	—	—
Radial Reload	Device	Reorder Code	Color Code	Staple Heights	Tissue Type	EES Product Code	EES Reload Color
		SIGRADMT	Purple	3 mm, 3.5 mm, 4 mm	Medium/Thick	CS40B/CS40G	Blue/Green
		SIGRADXT	Black	4 mm, 4.5 mm, 5 mm	Extra Thick	—	—


### TA™ Stapler With DST Series™ Technology

MEDTRONIC								ETHICON					
Shaft Length	Device	Reorder Code	Color Code	Reloads Available	Staple Height	Tissue Type	Max. Firings	EES Shaft Length	EES Product Code	EES Reload Code			
30 mm		TA30V3S	White	TA30V3L	2.5 mm	Vascular	8	30 mm	TX30V	XR30V			
		TA3035S	Blue	TA3035L	3.5 mm	Medium			TX30B	XR30B			
		TA3048S	Green	TA3048L	4.8 mm	Thick			TX30G	XR30G			
45 mm		TA4535S	Blue	TA4535L	3.5 mm	Medium	60 mm	TX60B	XR60B				
		TA4548S	Green	TA4548L	4.8 mm	Thick							
60 mm		TA6035S	Blue	TA6035L	3.5 mm	Medium				TX60G	XR60G		
		TA6048S	Green	TA6048L	4.8 mm	Thick							
90 mm		TA9035S	Blue	TA9035L	3.5 mm	Medium				90 mm	TL90	TR90	
		TA9048S	Green	TA9048L	4.8 mm	Thick							TLH90

### GIA™ Stapler With Tri-Staple™ Technology

MEDTRONIC								ETHICON		
Shaft Length	Device	Reorder Code	Color Code	Reloads Available	Staple Height	Tissue Type	Max. Firings	EES Shaft Length	EES Product Code	EES Reload Code
60 mm		GIA60MTS	Purple	GIA60MTC	3.0, 3.5, 4.0 mm	Medium/Thick	8	55 mm	TLC55 /NTLC55	TCR55
		GIA60XTS	Black	GIA60MTC	4.0, 4.5, 5.0 mm	Thick/Extra Thick			TCT55	TRT55
80 mm		GIA80MTS	Purple	GIA80MTC	3.0, 3.5, 4.0 mm	Medium/Thick		75 mm	TLC75 /NTLC75	TCR75
		GIA80XTS	Black	GIA80MTC	4.0, 4.5, 5.0 mm	Thick/Extra Thick			TCT75	TRT75

### GIA™ Stapler with DST Series™ Technology

MEDTRONIC								ETHICON		
Shaft Length	Device	Reorder Code	Color Code	Reloads Available	Staple Height	Tissue Type	Max. Firings	EES Shaft Length	EES Product Code	EES Reload Code
60 mm		GIA6025S	White	GIA6025L	2.5 mm	Vascular	8	55 mm	TVC55	TVR55
		GIA6038S	Blue	GIA6038L	3.8 mm	Medium			TLC55	TCR55
		SGIA6038S	Blue	GIA6038L	3.8 mm	Medium				
80 mm		GIA6048S	Green	GIA6048L	4.8 mm	Thick	75 mm	TCT55	TRT55	
		GIA8038S	Blue	GIA8038L	3.8 mm	Medium		TLC75	TCR75	
		GIA8038S	Blue	GIA8038L	3.8 mm	Medium				TCD75
		GIA8048S	Green	GIA8048L	4.8 mm	Thick		TCT75	TRT75	
100 mm		GIA10038S	Blue	GIA10038L	3.8 mm	Medium	100 mm	TLC10	TCR10	
		GIA10048S	Green	GIA10048L	4.8 mm	Thick		TCT10	TRT10	


†While offering a slimmer fixed anvil, reloads is a vascular reload that does not come with Tri-Staple™ technology.

© 2021 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.  
 †\*\* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.  
 07/2021-US-ST-2100120-[WF#2835774]

# STAPLING SOLUTIONS

## Product Reference


### EEA™ Circular Stapler with Tri-Staple™ Technology

MEDTRONIC						ETHICON	
Lumen Size	Device	Reorder Code	Color Code	Staple Height	Tissue Type	EES Lumen Size	EES Product Code
28 mm		TRIEEA28MT	Purple	3.0, 3.5, 4.0 mm	Medium/Thick	29 mm	CDH29A / CDH29P
		TRIEEA28XT	Black	4.0, 4.5, 5.0 mm	Thick/Extra Thick	29 mm	CDH29A / CDH29P
		TRIEEAXL28MT	Purple	3.0, 3.5, 4.0 mm	Medium/Thick	29 mm	ECS29A / ECS29B
		TRIEEAXL28XT	Black	4.0, 4.5, 5.0 mm	Thick/Extra Thick	29 mm	ECS29A / ECS29B
31 mm		TRIEEA31MT	Purple	3.0, 3.5, 4.0 mm	Medium/Thick	29 mm	CDH29A / CDH29P
		TRIEEA31XT	Black	4.0, 4.5, 5.0 mm	Thick/Extra Thick	29 mm	CDH29A / CDH29P
		TRIEEAXL31MT	Purple	3.0, 3.5, 4.0 mm	Medium/Thick	29 mm	ECS29A / ECS29B
		TRIEEAXL31XT	Black	4.0, 4.5, 5.0 mm	Thick/Extra Thick	29 mm	ECS29A / ECS29B

### EEA™ Circular Stapler with DST Series™ Technology

MEDTRONIC						ETHICON	
Lumen Size	Device	Reorder Code	Color Code	Staple Height	Tissue Type	EES Lumen Size	EES Product Code
21 mm		EEA2135	Aqua	3.5 mm	Medium	21 mm	CDH21A
		EEA21	Aqua	4.8 mm	Thick		
		EEAXL2135	Aqua	3.5 mm	Medium	21 mm	ECS21A
		EEAXL21	Aqua	4.8 mm	Thick		
25 mm		EEA2535	White	3.5 mm	Medium	25 mm	CDH25A
		EEA25	White	4.8 mm	Thick		
		EEAXL2535	White	3.5 mm	Medium	25 mm	ECS25A
		EEAXL25	White	4.8 mm	Thick		
28 mm		EEA2835	Blue	3.5 mm	Medium	29 mm	CDH29A
		EEA28	Blue	4.8 mm	Thick		
		EEAXL2835	Blue	3.5 mm	Medium	29 mm	CDH29 ECS29
		EEAXL28	Blue	4.8 mm	Thick		
31 mm	EEA31	Green	4.8 mm	Thick	29 mm	CDH29 ECS29	
	EEAXL31	Green	4.8 mm	Thick			
33 mm	EEA33	Yellow	4.8 mm	Thick	33 mm	CDH33A ECS33A	
	EEAXL33	Yellow	4.8 mm	Thick			

### EEA™ Hemorrhoid Stapler

MEDTRONIC						ETHICON	
Lumen Size	Device	Reorder Code	Color Code	Staple Height	Tissue Type	EES Product Code	
33 mm		HEM3335	Blue	3.5 mm	Medium	PPH03	
		HEM3348	Green	4.8 mm	Thick		